

Description

A food for special medical purposes.

Basecal® is a powdered protein free feed, containing carbohydrate, fat, vitamins, minerals, trace elements and the long chain polyunsaturated fatty acids (LCPs); Arachidonic acid (AA) and Docosahexaenoic acid (DHA).

Indications

Basecal is for use in the dietary management of inborn errors of protein metabolism requiring a low protein diet from 3 years of age.

Dosage and Administration

To be determined by the clinician or dietitian and is dependent on the age, body weight and medical condition of the patient.

> Basecal100 (21.5g) = 100kcal. Basecal200 (43g) = 200kcal.

Preparation Guidelines

For the standard dilution of 1kcal/ml:



Empty the contents of the sachet into the beaker.



Add cold water: For Basecal100 fill to a final volume of 100ml. For Basecal200 fill to a final volume of 200ml.



Replace lid and shake well for 10 seconds.



Basecal is now ready to use.

For best results, Basecal should be freshly prepared and where possible used immediately. Any unused Basecal should be refrigerated. Discard any remaining product 24 hours after preparation.

Tube feeding

Basecal can be used as part of a modular feed. Use as directed by your clinician or dietitian.

IMPORTANT NOTICE

Nutritional Information

		Per 100g	Per Basecal100 21.5g	Per Basecal200 43g
Energy	kJ	1952	420	839
	kcal	465	100	200
Fat	g	20.7	4.5	8.9
of which saturates	g	2.1	0.45	0.90
mono-unsaturates	g	13	2.8	5.6
polyunsaturates	g	4.7	1.0	2.0
Carbohydrate	g	69.8	15	30
of which sugars	g	9.5	2.1	4.1
Protein	g	0	0	0
Salt	g	0.49	0.11	0.21
Vitamins and min				
Vitamin A (RE)	μg	302	65	130
Vitamin D	μg	7.9	1.7	3.4
Vitamin E (aTE)	mg	9.3	2.0	4.0
Vitamin K	μg	26.5	5.7	11.4
Vitamin C	mg	69.8	15	30
Thiamin	mg	0.51	0.11	0.22
Riboflavin	mg	0.72	0.16	0.31
Vitamin B6	mg	0.60	0.13	0.26
Niacin mg (mg	g NE)	5.1 (5.1)	1.1 (1.1)	2.2(2.2)
Folic acid	μg	86	18.5	37
Vitamin B12	μg	1.3	0.28	0.56
Pantothenic Acid	mg	2.6	0.55	1.1
Biotin	μg	16.3	3.5	7.0

Nutritional Information

		Per 100g	Per Basecal100 21.5g	Per Basecal200 43g
Sodium	mg	195	42	84
	mmol	8.4	1.8	3.6
Chloride	mg	260	56	112
	mmol	7.3	1.6	3.1
Potassium	mg	465	100	200
	mmol	11.6	2.5	5.0
Calcium	mg	270	58	116
	mmol	6.8	1.5	2.9
Phosphorus	mg	195	42	84
	mmol	6.2	1.3	2.7
Magnesium	mg	47.4	10.2	20.4
	mmol	1.9	0.4	0.8
Iron	mg	4.7	1.0	2.0
Zinc	mg	4.7	1.0	2.0
Copper	mg	0.33	0.07	0.14
lodine	μg	79.1	17	34
Selenium	μg	14	3.0	6.0
Manganese	mg	0.47	0.10	0.20
Chromium	μg	13	2.8	5.6
Molybdenum	μg	22.3	4.8	9.6
Others				
Choline	mg	95.8	20.6	41.2
Inositiol	mg	51.2	11	22
DHA	mg	74	16	32
AA	mg	74	16	32

Allergen Declaration

Contains fish (tuna oil), milk (milk protein), soya (soya lecithin)

Storage

Store in a cool dry place. Sachets are designed for single use and should be used immediately once opened.

Pack Size/Weight

Basecal 100: $30 \times 21.5 \text{g} \oplus \text{sachets} = 645 \text{g}$ Basecal 200: $30 \times 43 \text{g} \oplus \text{sachets} = 1.3 \text{kg}$

Flavour

Unflavoured



